ISTRUZIONI D'USO

TALCUM S.A.L.F.

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TALCUM S.A.L.F.

Sterile, pyrogen-free talc powder for pleurodesis

CE 0373

TYPE OF DEVICE

Medical Device of Class IIb.

QUALITATIVE AND QUANTITATIVE COMPOSITION

Talcum S.A.L.F. is talc powder, a magnesium hydrated phyllosilicate [Mg₃Si₂O₇(OH)₂], free of asbestos, in conformity with the monograph 0438 of the current European Pharmacopoeia. The powder has a controlled particle size, the average particle size is 26 µm, in order to avoid systemic diffusion.

Talcum S.A.L.F. is packaged in an amber glass bottle and is available in the following containers:

50 ml bottle containing 2 g, 3 g, 4 g or 5 g of talc in single pack and / or 4 bottle pack.

30 ml bottle containing 1 g, 2 g, 3 g and talc in a single pack and / or 5 bottle pack.

INTENDED USE

Talcum S.A.L.F. is used for permanent pleurodesis. It is insoluble, sterile, non-pyrogenic powder, with controlled granulometry to avoid the risk of systemic spread through the pleura.

The pleurodesis is indicated in cases of massive or recurring effusions in benign and/or malignant pathologies or in the case of spontaneous pneumothorax, in order to prevent the accumulation of recurrent pleural fluid or air.

MECHANISM OF ACTION

Permanent pleurodesis intervention consists of administering into the pleural space, through a cannula or a catheter, a sclerosing substance, such as talc, to create a permanent adhesion of the two pleurae.

In physiological conditions, the pleural cavity is an almost virtual space delimited between two membranes, visceral and parietal, that slide on one another during breathing. The space is occupied by a small amount of liquid, that lubricates the membranes and facilitates the flow.

At the light microscopy the pleurae appear characterized by a single layer of mesothelial cells and connective tissue rich in collagen and elastin. The liquid is continuously formed and reabsorbed.

On the other hand, in pathological conditions there is the formation of either a transudation collection or a real pleuritic reaction intended as an inflammatory state of the serosa. Regardless of the etiology (malignant pleural effusions with different primitivities, benign idiopathic pleural effusions, secondary benign pleural effusions), a chronic and excessive accumulation of liquid occurs in the cavity, which compromises the lung function.

The therapeutic approach consists in removing the accumulated fluid and introducing talc through a catheter to generate a permanent adhesion of the two pleurae.

The mechanism of action of talc in the pleural cavity is the direct consequence of its chemical-physical properties.

Talc is insoluble powder, that, if administered thoracoscopically, fills the pleural space of the patient, thus favouring the adhesion of the two membranes. Adhesion is favoured by the controlled particle size; they remain permanently occupying the cavity and do not diffuse through the pleura due to their size (systemic diffusion). From a chemical point of view, talc is a magnesium hydrated phyllosilicate consisting of tetrahedral units repeated in layers on a regular basis, that form OH-Si bonds between them.

In pathological conditions the pH inside the pleural cavity becomes acidic and this determines a change in the structure of the molecule; in an acidic environment the hydroxyl group (OH) reacts with the trivalent iron present in the pleural space, determining oxidative stress and adhesion of the membranes.

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The insufflation of talc occurs after full aspiration of the pleural liquid and removal of any adhesions, through a video-guided administration. The amount of talc to be used is established by the doctor in relation to the diagnosis, the patient's age and health conditions. The dosage varies from 2 to 6 grams for pleural effusion and generally it does not exceed 2 grams for spontaneous pneumothorax.

Talcum S.A.L.F. can be introduced in the pleural cavity between the two membranes by pulverization (poudrage) thoracoscopically, or the powder can be mixed in advance with a sodium chloride saline solution 0.9% and then injected into the pleural cavity (slurry) through a drainage tube or through a thoracic tube. In this case, every gram of powder is suspended with 10 ml of sterile sodium chloride solution 0.9%.

Regardless of the technique chosen, the procedure requires that the pleural space is free of liquids. Therefore, the pleural fluid present is drained before proceeding.

To use the medical device, please, follow the approved operating protocol.

Below are some operating methods:

PULVERIZATION (POUDRAGE)

- Completely remove the aluminum capsule from the bottle
- Connect the bottle to a nebulizer
- Through a cannula, spread talc homogeneously inside the pleural space to obtain a "nebulization" of the product
- Do not apply an excessive amount of talc
- To obtain a homogeneous distribution of talc in the pleural cavity, repeat the operation modifying the orientation of the cannula

- During the application, gently swing the bottle to prevent talc from adhering to the internal walls
- Talc should not be sprayed only at one application site
- To correctly apply 2 g of talc, it takes about 2 minutes with a maximum frequency equal to 1 nebulization per second. It is normal that at the end of the operation the bottle still contains a residue of talc without this compromising the effectiveness of the application.

SUSPENSION (SLURRY)

Talcum S.A.L.F. can also be applied after suspension. The use of 10 ml of sterile solution of sodium chloride 0.9% for every gram of talc is recommended. The suspension has to be used immediately after its preparation. Thoroughly shake the bottle to avoid agglutination. Apply with a drainage tube or a washing syringe of 50 ml and inject into the thoracic cavity through a chest tube.

It is recommended to change regularly the position of the patient, in order to favour a homogeneous distribution of the suspension.

CONTRAINDICATIONS

Talcum S.A.L.F. must not be used:

- In patients who cannot be subjected to thoracoscopy;
- In patients with mechanical abrasion of the pleura, in order to avoid systemic migration;
- In case of pregnancy and during breastfeeding.

INTERACTIONS

Talc is incompatible with quaternary ammonium compounds.

ADVERSE REACTIONS

In consequence of pleurodesis with talc, the following undesirable effects may occur:

- Fever
- Pleural or chest pain
- Breathing difficulties (dyspnoea)

Clinical experience and literature highlight the effect of the different particle size in relation to the undesirable effects described: smaller particles migrate more easily. In particular, talc containing a quantity equal to 50% of particles smaller than 10 µm (mixed talc) causes greater reactions of systemic inflammation and hypoxemia (insufficient oxygen content in the blood) compared to talc containing less than 50% of particles smaller than 20 µm (graded talc). The particles smaller than 10 µm in sterile Talcum S.A.L.F. are less than 30%.

The influence of this dimensional factor is important in relation to the decrease in the appearance of acute respiratory distress syndrome or other complications.

SPECIAL WARNINGS AND APPROPRIATE PRECAUTIONS FOR USE

- The device is to be used only by qualified healthcare personnel working with aseptic technique in compliance with the relative operating protocol;
- Do not use, if the primary container is damaged;
- Use immediately after opening. The device is intended for a single use only, any residue cannot be used and should be eliminated due to the risk of contamination, which causes loss of sterility;
- The device is not to be administered orally;
- Following pleural powdering, a thoracic X-ray could highlight mistakenly the talc remained in the pleural cavity as a nodule, cancer or cyst;
- The use of FANS following chemical pleurodesis is controversial and could reduce the effectiveness of the technique: their use should be carefully evaluated by the specialist doctor who is treating the patient.

WARNINGS AND PRECAUTIONS IN RELATION TO THE VALIDITY OF THE DEVICE

Check the expiry date indicated on the primary container of the device. The expiry date refers to the product in intact packaging.

ATTENTION: do not use the device after the indicated expiry date.

SHELF LIFE

36 months referred to the product in unopened packaging, correctly stored.

After the first opening: the device should be used immediately.

SPECIAL PRECAUTIONS FOR STORAGE

Store in the original package, away from heat sources. No special storage conditions with respect to temperature are required. Do not refrigerate or freeze. The sterility of the product is guaranteed only when unopened. **KEEP OUT OF THE REACH AND SIGHT OF CHILDREN.**

NATURE AND CONTENTS OF THE PACK

The medical device is packaged in an amber glass bottle. It is available in a single and/or multiple package of 4 or 5 bottles.

SPECIAL PRECAUTIONS FOR DISPOSAL AND HANDLING OF THE DEVICE

Avoid release to the environment. Any unused product or waste material should be disposed of in accordance with local requirements.

MANUFACTURER

S.A.L.F. S.p.A. Laboratorio Farmacologico - Via Marconi, 2 - Cenate Sotto (BG), Italy

SYMBOLS

The symbols stated on the container and relative descriptions are listed below:

	Medical Device Code
	LOT Number
	Expiry Date
	Protect from direct sun light
	Do not reuse / Single use only. Reuse can lead to the risk of contamination and loss of sterility
	Do not re-sterilise
	Pay attention to the instructions for use!
	Sterilized using irradiation
	Do not use if the immediate container is damaged
	Manufacturer

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Talcum S.A.L.F. é embalado num frasco de vidro âmbar e está disponível nos seguintes formatos: